NORTHERN LIGHT TECHNOLOGIES

Website: www.NorthernLightTechnologies.com
Tel: (514) 335-1763 Toll Free: 1-800-263-0066
Email: info@NorthernLightTechnologies.com

Address: 8971 Blvd Henri-Bourassa W., Montreal, QC, H4S 1P7

INSTRUCTIONS FOR USE

The **TRAVelite II** is designed to bring spring and summer outdoor light levels into your home or office for years to come. Use the lamp for **30 minutes each day**, preferably first thing in the morning. Sit in front of the lamp, **with your eyes 12 inches from the center** of the diffuser. This will ensure that you are exposed to 10,000 Lux. Do not look directly at the light, but you may glance at it periodically. Install stand as shown on sketch. You can use the unit either in a vertical or horizontal position (see above). Keeping regular sleep/wake-hours, and exercising (outdoor when possible) are also very beneficial.

PROBLEMS

1) UNIT DOES NOT LIGHT

If the unit does not light after the switch is turned ON, please check that the plug is well inserted into the outlet. If it still fails to turn on, try other outlets. Should the problem persist, remove the plastic diffuser by gently unscrewing the four nylon screws. Remove the bulbs and reinsert them ensuring that both tubes are properly seated in the sockets.

2) OTHER PROBLEMS

Bright light can cause eye and/or skin irritation and headaches. This is most prevalent in people using

light sensitizing medications, or people with blue or green eyes, blonde or red hair, fair or sensitive skin. These problems are usually resolved by gradually increasing use of the lamp. At the first sign of skin or eye irritation, stop using the lamp and wait until the irritation disappears.

Problems caused by bright light can usually be overcome by first sitting 2 feet away for 2 minutes and gradually increasing the light over two weeks, until you are 12 inches for at least 20 minutes.

If you still experience problems, call us. We will try to help you.

WARNING: *If you have a health problem, if you have been told you must wear sunglasses in sunlight, or if you have a history of eye problems, please consult a health-care professional first. MEDICATIONS such as lithium, melatonin, tetracycline, St. John's Wort, acne creams, creams with retinoic acid and some chemotherapy drugs can make you very sensitive to light.*

MAINTENANCE

The lamp uses two 18 watt LED tubes with a 30,000 hour duration. The tubes must be purchased from Northern Light Technologies. To replace the tubes, gently remove the 4 plastic diffuser screws, put in new tubes and put the diffuser back in place.

WARRANTY

Northern Light Technologies will repair or replace a defective lamp at no charge for a period of 7 years. Shipping costs to our plant must be paid by you. You **MUST** call us at **1-800-263-0066** before returning a unit. Warranty does not cover the tubes.

GUARANTEE

For a period of 60 days after purchase we will refund the total paid, less \$35 for shipping and restocking. If the lamp was not purchased from us, please see vendor for return policy.

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